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e-mail: info@eurosets.com - homepage: www.eurosets.com

WaterLily™ Wound Therapy

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**WaterLily Vacuum Generator TS
for Negative Pressure Wound Therapy**
Multiuse vacuum generator for negative pressure wound therapy

CE₀₁₂₃

USER INSTRUCTION MANUAL

Index

General warnings.....	59
Product description.....	60
Intended use.....	61
Contraindications.....	63
Precautions.....	63
Warnings	64
Ambient storage and transport conditions	65
Ambient conditions of use	65
Technical description of the device.....	65
Electrical/Electronic	65
Battery life.....	65
Pneumatic	65
Suction pressure.....	65
Suction flow.....	66
Power source management.....	66
Recharging the battery.....	66
Rechargeable internal battery low warning.....	67
Rechargeable internal battery exhausted warning.....	67
Cleaning.....	68
Periodic maintenance	68
Internal rechargeable battery life cycle	68
Other malfunctions that may occur during normal device operation.....	68
Periodic maintenance schedule	68
Device lifetime	68
Preparing for use	69
WaterLily TS Canister installation procedure.....	69
Vacuum generator attachment, support and carrying systems	69
Switching on	71
Checking the alarms.....	71
ALARM CHECK: functional test of the visual and acoustic alarm system	71
Activating NEW THERAPY.....	72
COMFORT and RAPID mode.....	73
Symbols used on the operating page.....	75
Basic programming of the suction pressure	76

Replacing the WaterLily TS Canister	76
Switching off the device	76
Switching on the device and activating CONTINUE THERAPY	76
Using the advanced functions	77
 Changing the parameters in the main menu	78
1 - ASPIRATION.....	78
Intermittent.....	79
Tube flushing.....	80
2 - MANOMETER	80
Scale	80
3 – TIMER	82
Activating the timer and setting a reminder.....	82
4 – ALARMS	83
Automatic alarm reset	83
Volume	83
5 – HISTORY.....	84
EVENT SUMMARY	85
EVENT DETAILS.....	85
EVENT DOWNLOAD.....	86
6 – KEY LOCK	88
KEY LOCK	89
7 - SETTINGS	90
7.1 - LANGUAGES	91
7.2 - DATE / TIME	91
7.3 - DISPLAY	92
Keyboard beep	92
Screensaver.....	93
7.4 - ALARM CHECK.....	93
7.5 - RESET.....	94
7.6 - MAINTENANCE.....	95
7.7 - UPDATE	95
 Alarms	96
Information signals	102
 Alarm and information signal priorities in case of overlapping.....	103
 Device classification according to the standards and relative limits of use	104
Electromagnetic compatibility.....	105
ELECTROMAGNETIC EMISSIONS	105
ELECTROMAGNETIC IMMUNITY.....	106
RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND WATERLILY VACUUM GENERATOR TS	108
Warranty conditions	108
Explanation of the symbols used on the device labels.....	109

General warnings

- *Carefully read the instructions before use.*
- *The device may only be applied by qualified and suitably trained healthcare workers.*
- *This device is composed of precision components. Do not drop it or subject it to strong impact. Never attempt to disassemble it.*
- *No device hardware and/or software modification is permitted unless performed by the manufacturer.*
- *The device must be used in accordance with these instructions. EUROSETS cannot be held responsible for damage deriving from improper use.*
- *To prevent the risk of electric shock or fire or device malfunctioning:*
 - *Do not wet the device.*
 - *Do not use the device in the shower.*
 - *Position the device in such a way that it cannot fall into water or other liquids (e.g. basin). Should it accidentally fall into water, do not touch the device and immediately pull the plug from the power socket.*
 - *Do not expose the device to the steam generated by kettles, atomizers, irons, etc.*
 - *Do not expose the device to direct heat sources (radiators, fireplaces, etc.).*
 - *Only use the power supply/battery charger provided.*
- *Do not leave the device unattended (including the power supply/battery charger) where there are children, persons with reduced cognitive abilities and animals.*
- *To isolate the device from the main electric line, disconnect the power supply from the mains.*
- *Do not introduce the device into a hyperbaric chamber as it is not designed for operation in an oxygen-rich environment.*
- *Electromagnetic interference:*
 - *Do not use the device in the vicinity of equipment not sufficiently screened from electromagnetic emissions higher than those indicated in IEC EN 60601-1 (2007).*
 - *Reciprocal interference may occur between this device and wireless communication devices, such as domestic wireless devices, mobile phones, cordless phones and their base radio stations, walkie-talkies. In the event of reciprocal interference between devices, move the communication devices away from this device. The safety distance is indicated in the paragraph DEVICE CLASSIFICATION ACCORDING TO THE STANDARDS AND RELATIVE LIMITS OF USE.*
- *Reusable device. Do not discard after use.*
-  *Do not dispose of WEEE (waste electrical and electronic equipment) as solid household waste but collect it separately; the device can be returned to the distributor/manufacturer.*
The hazardous substances contained in electrical equipment may be carcinogenic to man if dispersed into the environment.

Product description

WaterLily Vacuum Generator TS is a multiuse device for negative pressure wound therapy and constitutes the medical device of the WaterLily System.

The WaterLily System comprises:

- Multiuse **WaterLily Vacuum Generator TS**.
- Single-use **WaterLily TS Canister** for fluid collection.
- Tube with suction head or, alternatively, suction tube connected to a Jackson-Pratt drain, both single-use.
- Single-use gauze or porous foam wound dressing.
- Single-use adhesive drape.

The **WaterLily Vacuum Generator TS** must only be used connected to the single-use **WaterLily TS Canister** and in conjunction with the suction tube, wound dressing and adhesive drape.

Carefully follow the instructions for use provided with the canister, suction tube, wound dressing and adhesive drape before starting to use the **WaterLily Vacuum Generator TS**.

The single-use **WaterLily TS Canister** is equipped with a hydrophobic sterilizing membrane that acts as overflow protection and prevents contamination of the **WaterLily Vacuum Generator TS**.

Caution

*The moment the overflow protection system is activated, the vacuum generator stops supplying vacuum useful for the negative pressure wound therapy.
Replace the canister to resume the therapy.*

Illustration of the WaterLily system with indication of its main components:

- 1) Single-use collection canister
- 2) LCD display (touchscreen)
- 3) Multiuse vacuum generator
- 4) USB port
- 5) Green and yellow/red LED
- 6) On/off button
- 7) Power supply connection port
- 8) Tube with suction head that is connected to the patient's wound dressing

The following devices are an integral part of **WaterLily Vacuum Generator TS**:

9-10) Two modular straps adjustable in length to secure the device to the edge of the bed, carry it by hand or slung over the chest or shoulder

11) USB pen drive (type: "Eurosets")

12) Power supply/battery charger (type: FW7362M/18)

Intended use

WaterLily Vacuum Generator TS connected to the single-use **WaterLily TS Canister** and in conjunction with the suction tube, wound dressing and adhesive drape is suitable for performing negative pressure wound therapy in a hospital or domestic environment.

WaterLily Vacuum Generator TS is suitable for application of suction (negative pressure) to the patient's wound in order to:

- increase local blood irrigation;
- reduce oedema;
- stimulate the formation of granulation tissue;
- stimulate cell proliferation;
- eliminate soluble healing inhibitor agents from the wound;
- reduce the bacterial charge;
- bring together the wound edges.

WaterLily Vacuum Generator TS is suitable for treatment of the following wounds:

- chronic ulcers on the lower limbs;
- diabetic foot ulcers;
- open abdomen including explored fistulas;
- pressure ulcers;
- skin graft take;
- infected sternal wounds;
- surgical dehiscence;
- traumatic wounds.

The recommended maximum negative pressure (suction) varies between 80 and 120 mmHg depending on the therapeutic protocol at the discretion of the attending physician.

The negative pressure (suction) generated by **WaterLily Vacuum Generator TS** can be set by the operator based on the therapeutic protocols and is between **20** and **200 mmHg** with increments of 5 mmHg.

The negative pressure can be applied continuously or intermittently with time control according to the therapeutic protocol.

WaterLily Vacuum Generator TS has been designed for basic use or for advanced use by accessing the **MENU** functions and changing the preset basic parameters (language, suction pressure unit of measure, etc.).

WaterLily Vacuum Generator TS has a data logger function, i.e. you can record the set parameters and the operating events of the device, which can be:

- viewed on the display;
- downloaded to the USB pen drive provided in CSV format and viewed on a PC with Microsoft Excel.

Caution

Any negative pressure wound therapy application outside the scope of application specified by the manufacturer is at the discretion and under sole responsibility of the attending physician.

Contraindications

- Do not apply the wound dressing directly on exposed organs, blood vessels, nerves, ligaments and tendons.
- Do not use the negative pressure wound therapy:
 - on an exposed anastomotic site,
 - if you have even the slightest suspicion that the wound may be malignant,
 - in proximity of the cranial nerve in order to minimise the risk of bradycardia,
 - in proximity of a spinal cord injury, which might cause undesired stimulation of the sympathetic nervous system with consequent autonomic dysreflexia,
 - in the presence of necrotic tissue with a eschar on the wound,
 - in the presence of osteomyelitis,
 - in the presence of unexplored fistulas and enteric fistulas,
 - in the presence 3rd degree burns with a eschar on the wound.

Precautions

- Before application of the negative pressure wound pressure therapy, carefully assess the weight and size of the patient. For patients at risk of fluid depletion or dehydration, it is suggested to set a low initial negative pressure and check the volume exuded.
- Particular attention must be paid in using the hypobaric therapy on the following patients:
 - patients with coagulation disorders,
 - patients that have taken or have recently taken anticoagulants and/or antiaggregants,
 - patients that have recently suffered acute or surgical traumas,
 - patients with weak or infected blood vessels,
 - patients with infected wounds,
 - patients with sharp parts in the wound.

The risk of haemorrhage increases considerably in such patients and it is advisable to take the following precautions:

- in case of surgery, postpone application of the negative pressure wound therapy for at least 24 hours,
- when applying the negative pressure wound therapy, frequently check the fluid collection canister, the suction tubes and the drains for any signs of bleeding,
- be extremely careful when changing the dressing, as this is when the risk of haemorrhage is the highest.

Warnings

Warnings

The absence of vacuum in the dressing might determine build-up of exudates triggering infectious and inflammatory processes. If the absence of vacuum extends for several days, necrotic eschars may form. If suction is interrupted for more than 2 hours, remove the dressing after irrigating it with saline solution.

Use in a hyperbaric chamber

For patients connected to the WaterLily system who require hyperbaric treatment, follow the instructions below:

- 4. Do not introduce the vacuum generator into a hyperbaric chamber as this device is not suitable for use in an oxygen-rich environment.*
- 5. Before the patient enters the hyperbaric chamber, switch off the vacuum generator and disconnect the single-use canister (see the procedure in the paragraph "Replacing the WaterLily TS Canister").*
- 6. When the patient exits the hyperbaric chamber, reconnect the canister to the vacuum generator and resume the therapy using the CONTINUE THERAPY function (see the procedure in the paragraph "Turning on the device and activating CONTINUE THERAPY").*

Magnetic resonance imaging (MRI) environment

For patients connected to the WaterLily system who require magnetic resonance imaging (MRI), follow the instructions below:

- 4. Do not introduce the vacuum generator into a magnetic resonance imaging environment as this device is not suitable for use in such an environment.*
- 5. Before the patient enters the magnetic resonance imaging environment, switch off the vacuum generator and disconnect the single-use canister (see the procedure in the paragraph "Replacing the WaterLily TS Canister").*
- 6. When the patient exits the magnetic resonance imaging environment, reconnect the canister to the vacuum generator and resume the therapy using the CONTINUE THERAPY function (see the procedure in the paragraph "Turning on the device and activating CONTINUE THERAPY").*

Use in case of defibrillation

For patients connected to the WaterLily system who require defibrillation in the area involving the dressing, remove the dressing before performing the treatment.

Ambient storage and transport conditions

- Store the device in a place where the temperature is between 0° and 40°C; if exposed to low temperatures, the battery may spontaneously discharge.
- Store the device in a dry place not subject to humidity (relative humidity <90%).
- Store and transport the device in the box in which it was delivered.

Ambient conditions of use

- Use the device in a place where the temperature is between 5° and 35°C, relative humidity < 90% and with an atmospheric pressure between 70 KPa and 106 KPa.

Technical description of the device

Electrical/Electronic

- Device with circuit board and microprocessor
- 4.3" backlit colour LCD-TFT display with resistive touchscreen
- Power supply with 2 alternative sources:
 - Internal lithium ion battery pack (Li-ion Polymer Cell) (4 x 3.7V), rechargeable only with the power supply/battery charger provided,
 - Power supply/battery charger provided:
Input: 100-240VAC / 50-60 Hz / 700 mA
Output: 18VDC / 1.66 A
- The priority of use of the energy sources is as follows:
 - 1 - Mains
 - 2 - Rechargeable battery

Battery life

During typical use, the rechargeable battery can supply energy for up to 14 hours of use before having to be recharged.

Battery life depends on many variables, such as:

- *conditions of the wound and the dressing, for example, seal of the adhesive film on the skin,*
 - *number of battery discharge and charge cycles,*
 - *battery discharge and charge conditions,*
 - *device storage conditions,*
 - *date of manufacture of the device.*
- The battery life indicated is hence purely approximate.*

Pneumatic

Suction pressure

- Suction pressures between **20** and **200 mmHg**.
- Suction pressure unit of measure: **mmHg**.
- Suction pressure unit of measure selectable in **ADVANCED USE: cm H₂O, mBar, Pa x 100**.
- The suction pressure is read and adjusted by the electronic sensor fitted in the device with tolerances of **-5 /+15 mmHg**.

Suction flow

- The device is equipped with an automatic suction flow rate adjustment system.
- The automatic adjustment system is controlled by the software installed in the device and allows obtaining the ideal suction flow for the wound seal conditions and the set suction pressure.


Power source management

Recharging the battery

- ***The internal rechargeable battery is supplied exhausted.***
- ***If the device needs to be used not connected to the mains, check the battery charge before starting the therapy.***

- When the rechargeable battery is flat, the device can be powered directly from the mains by means of the power supply/battery charger which at the same time as supplying power also starts a new recharging cycle.
- The device can be connected to the power supply/battery charger both when it is on and operational and when it is off and not operational.
- The rechargeable battery must be recharged using the power supply/battery charger provided.
- Connect the battery charger to the vacuum generator port (right-hand side) and then to the mains at a point easily accessible to the operator.
- While the battery is charging, the following is shown on the display:

Device on and operational with the power supply/battery charger operational (the battery charge level varies dynamically).

Device off and not operational with the power supply/battery charger operational.

- It takes 3-4 hours at a temperature of 20°C + 5°C to recharge the lithium ion battery pack.
- When battery charging is complete, the following is shown on the display:

Device on and operational with the power supply/battery charger operational with power supply function only.

Device off and not operational with the power supply/battery charger no longer operational.

Warnings

- *There is a risk of tripping over the power supply/battery charger cable. Make sure that the cable is away from areas where people pass through.*
- *The length of the power supply/battery charger cable is such that it may cause strangulation/suffocation.*
- *The battery cannot be accessed or repaired by the user.*

Rechargeable internal battery low warning

The battery charge is shown at the top right of the screen.

The battery low symbol appears when the battery charge is at about 1/5 of its total capacity.

When the battery is low, a battery charge warning and a yellow blinking triangle appear on the screen accompanied by an acoustic signal and the yellow LED blinking.

Rechargeable internal battery exhausted warning

The battery exhausted symbol appears when the battery has completely discharged.

The red triangle indicates that the battery is exhausted.

When the battery is exhausted, a battery charge warning and a red blinking triangle appear on the screen accompanied by an acoustic signal and the red LED blinking.

Caution

- *To maintain battery life and perfect efficiency, it is advisable to recharge it when the device signals that the battery is low or exhausted.*
- *If the display does not turn on during the first few minutes of recharging, it does not mean that there is a malfunction, but only that the battery is exhausted.*
- Leave the device connected to the power supply/battery charger connected to the mains and the display will come on after a few minutes (it takes 3-4 hours to fully recharge the battery).*

Cleaning

- The outer casing is to be cleaned/disinfected by healthcare workers using a cloth soaked in a preferably water-based, wide-spectrum disinfectant solution (e.g. chlorhexidine), wringing it out well.

Warning

- Do not use chlorine-based disinfectant solutions as they will oxidise the metal parts.*
- Alcohol-based disinfectant solutions should only be used occasionally (once a month) and only with concentrations in water of not more than 70% in order to prevent damaging the plastic parts of the device.*
- The cleaning/disinfection operations should always be carried out before transferring the device to another patient and may be carried out during treatment cycles on the same patient.*
- Do not let liquids drip into the device during the cleaning operations.*
- Do not spray or pour disinfectant directly onto the device.*

Periodic maintenance

Internal rechargeable battery life cycle

The internal rechargeable battery has been tested to withstand at least 500 charge/discharge cycles in the intended conditions of use (see the paragraph **Ambient conditions of use**).

Nonetheless, because of the composition of this electronic component, the number of charge/discharge cycles may vary up or down. In the event of poor or no battery life, contact Eurosets Technical Service or an authorised agent.

Other malfunctions that may occur during normal device operation

In the event of electrical, electronic, pneumatic or mechanical faults, or if you suspect that liquids or solids have been sucked into the vacuum pump, contact Eurosets Technical Service or an authorised agent.

Periodic maintenance schedule

Device functioning needs to be checked at least every 3 years (36 months); this procedure includes the replacement of the internal rechargeable batteries. The healthcare facility is responsible for requesting technical service by the manufacturer or an authorised agent.

Device lifetime

Provided that the above mentioned maintenance schedule and procedures are observed, the device lifetime is 10 years (120 months) from the date of delivery to the end customer/healthcare facility.

Caution

The essential safety requirements of the device are not guaranteed in case of:

- maintenance (including battery replacement) and hardware and software modifications not carried out by the manufacturer or an authorised agent;*
- failure to observe the periodic maintenance schedule and to request technical service by the manufacturer or an authorised agent;*
- using the device beyond its lifetime of 10 years (120 months) from the date of delivery to the end customer/healthcare facility.*

Preparing for use

WaterLily TS Canister installation procedure

After connecting the single-use canister to the wound suction tube/drain and after placing it outside the sterile field, install it on the **WaterLily Vacuum Generator TS**.

The canister has two slots (Fig. 1 and Fig 2) for horizontally securing it on the vacuum generator and a snap hook (Fig. 3 and Fig. 4) for hooking it onto the vacuum generator (Fig. 4).

Connect the canister to the vacuum generator by inserting the canister into the guides (1), move it into vertical position with the wall against the vacuum generator (2), and then lower the hook until you hear a click (3).


Fig. 1


Fig. 2


Fig. 3

Fig. 4

Vacuum generator attachment, support and carrying systems

Two modular straps are supplied together with the device and can be adjusted to the desired length. The short strap (Fig. 5) can be used to secure the device to the edge of the bed (Fig. 6) or to carry it by hand (Fig. 7).

Fig. 5


Fig. 6

Fig. 7

The long strap (Fig. 8) can be used to carry the device slung over the chest (Fig. 9) or shoulder (Fig. 10).

Fig. 8

Fig. 9

Fig. 10

Hook the chosen strap onto the metal rings on the upper part of the device (Fig. 11).


Fig. 11

A protective neoprene cover for carrying the device around is also provided. Insert the device in the cover and lower the Velcro flap to ensure proper closure.

Switching on

To switch on the vacuum generator, press and hold the  button for a few seconds.

If the device switches on properly, an acoustic signal (beep) will sound, the green LED  will come on, and the home page followed by the start page will appear:

The vacuum generator is equipped with a touchscreen through which the operator controls the device. The device is supplied with the following factory settings:

- Communication language : **English**
- Suction pressure unit of measure: **mmHg**
- Date/time: **Central European Time (CET)**.

To change the above settings, refer to paragraphs **7.1 – LANGUAGES**, **2 – MANOMETER** and **7.2 – DATE/TIME**, respectively.

Checking the alarms

ALARM CHECK: functional test of the visual and acoustic alarm system

The **ALARM CHECK** must be carried out by healthcare workers the first time the device is switched on and before transferring the device to another patient following the procedure below:


After turning on the device and the start page has appeared:

press the **MENU** button to access the main menu, press the **SETTINGS** pane to access the **SETTINGS** menu, then press the **ALARM CHECK** pane.

The following page will appear:


accompanied by emission of the sequence of acoustic alarms and the following LEDs coming on in sequence:
GREEN ⇒ RED ⇒ YELLOW.

During the test, check that:

- the speaker emits the acoustic alarms in sequence,
- the LEDs come on in the 3 different colours.

At the end of the test, the following page will appear:

indicating that the test is complete.


If the acoustic and visual alarms have been emitted correctly, press **OK** to return to the **SETTINGS** menu, press **ESC** to return to the main menu, and press **ESC** again to return to the start page.


If any errors occur, switch off the device and contact Eurosets Technical Service or an authorised agent.

Activating NEW THERAPY


Press **NEW THERAPY** on the start page

and the following page will appear:


Check that the date and time are correct. If not, press the **NO** button to return to the start page:


next press **MENU** then **SETTINGS** and finally **DATE/TIME**.

Select the up or down arrows to change the date and time, then press **OK** to confirm.

Once the date and time have been set, press **OK**. The **SETTINGS** menu will automatically appear. Press **ESC** to return to the start page.

Press **NEW THERAPY** on the start page:

then select the desired suction mode: **COMFORT** or **RAPID**.

COMFORT and RAPID mode

The vacuum generator can be set to two suction modes:

- **COMFORT**
- **RAPID**

These two modes differ between them in the suction flows applied to the wound.

COMFORT: medium suction flow for the first few minutes from activation with a gradual increase if the set pressure is not reached up to the maximum suction flow.

INDICATIONS: moderately-sized wounds; less painful to the patient.

RAPID: Maximum suction flow right from activation.


INDICATIONS: large wounds.

UNDESIRED EFFECTS: more painful to the patient.

Note


After the first few minutes of application, suction flow control is automatic and independent of the values initially set.

After selecting the desired mode, **COMFORT** or **RAPID**, the page for entering the patient's initials appears. After entering the patient's initials (identification), press **OK** to continue and start the therapy. Press **ESC** to cancel.


After confirming with **OK**, the operating page appears. The words **SEAL TEST** will automatically disappear once the set suction pressure has been reached.

A two-tone acoustic signal will also be sounded when the set suction pressure has been reached.


Symbols used on the operating page

1	Dynamic GEAR symbol: indicates when the vacuum generator is running (when the gears move, the generator is in suction)
2	Left arrow button : press this button to decrease the set vacuum
3	Calendar : indicates the current date and time
4	Set pressure and relative unit of measure : indicates the set pressure and the relative unit of measure
5	Right arrow button : press this button to increase the set vacuum
6	BATTERY and POWER PLUG symbol: indicates the charge status of the internal rechargeable battery and whether it is on charge (indicated by the power plug symbol)
7	Dynamic LEAK bar: indicates the entity of air leakage on the dressing
8	PADLOCK symbol: if shown, it indicates that the touchscreen buttons are locked
9	TORTOISE or HARE symbol: indicates whether the therapy is being performed in COMFORT mode (tortoise) or RAPID mode (hare)
10	MENU button: press and hold the button for a few seconds to access the main menu
11	SEAL TEST symbol: indicates that the device is in the process of reaching the negative pressure in the dressing and has not yet sealed the dressing
12	Dynamic vacuum bar : dynamic indication of the negative pressure on the dressing
13	Alarm present symbol: appears when an alarm is present and indicates its priority (in the example, a high priority alarm)
14	Alarm silenced symbol: appears when an alarm is silenced

Basic programming of the suction pressure

Repeatedly press  or  depending on whether you want to increase or decrease the suction pressure. Each time an arrow is pressed, the suction pressure decreases or increases by 5 mmHg. Once you have set the new pressure value, a page appears asking you to confirm:

Press **YES** to confirm or **NO** to cancel.

The last suction pressure set remains in memory even if the device is switched off (it switches on again in **CONTINUE THERAPY** mode) and/or the battery goes completely flat.

Replacing the WaterLily TS Canister

To replace the single-use canister when it is full, disconnect it from the wound suction tube/drain and remove it from the vacuum generator by raising the hook and partially turning it anticlockwise and then extracting it upwards.

Switching off the device

Press and hold down the  button until the device switches off.

While the device is switching off, you will hear an acoustic signal (beep) and the following will appear on the display:

Switching on the device and activating CONTINUE THERAPY

Note

*When activating the **CONTINUE THERAPY** function, the preset parameters will be maintained and chronological recording of the events continued.*


Press and hold the  button for a few seconds until the home page followed by the start page appears:

Press **CONTINUE THERAPY** on the start page

and the following page will appear:

On this page, press **COMFORT** or **RAPID** to continue the previously interrupted therapy.

After selecting the desired mode, the operating page appears where the previously set parameters and pressures will have remained unchanged.


Using the advanced functions

To switch on the vacuum generator, press and hold the  button for a few seconds.


If the device switches on properly, an acoustic signal (beep) will sound, the green LED will come on, and the home page followed by the start page will appear:

Press and hold the **MENU** button on the start page for a few seconds:

The main menu containing the device function options appears.

1	ASPIRATION: suction mode settings
2	MANOMETER: unit of measure settings
3	TIMER: visual and acoustic reminder settings
4	ALARMS: alarm settings
5	HISTORY: page to view and download the events register of the therapy in progress
6	KEY LOCK: page to lock/unlock the buttons
7	SETTINGS: page to set the device configuration

Press **ESC** to return to the start page.

Note

During navigation and programming, the device continues performing all its functions, displaying any alarms that are generated.


Changing the parameters in the main menu

1 - ASPIRATION

- Intermittent operation adjustment
- **TUBE FLUSHING** function

Press the **ASPIRATION** pane to access the page that includes intermittent operation adjustment and the **TUBE FLUSHING** function.

Intermittent

Allows setting timed intermittent operation, i.e. generating and maintaining the set suction pressure according to the times determined by the operator, for example:

ON TIME: 5 min

OFF TIME: 1 min

This means that a cycle is activated where:



- the device maintains the pressure set by the operator for 5 minutes;
- the device goes into standby for 1 minute and allows the suction pressure to drop to a minimum value of 30 mmHg;
- the 5-minute phase is resumed;
- the 1-minute phase is resumed;
- etc.

Caution

- *The maximum pressure value achievable during the ON time is that set in main programming.*
- *The minimum pressure value achievable during the OFF time will not be less than 30 mmHg.*
- *The ON and OFF values are settable with increments of 1 minute (minimum value settable 1 minute).*

To activate this function, press on the switch-shaped button until it turns green and the word **YES** appears. When the function is active, you can change the **ON** and **OFF** times by pressing on the relative arrows. When setting is complete, confirm with **OK**. The main menu will automatically appear.

Caution

- *The intermittent ON/OFF values set are shown on the display.*
- *ON or OFF is shown intermittently on the display according to the operating phase in which the device is.*

Tube flushing

Allows setting suction tube/drain flushing with air preventing possible formation of clots/occlusions in the tube. To activate this function, press on the switch-shaped button until it turns green and the word **YES** appears. When setting is complete, confirm with **OK**. The main menu will automatically appear.

2 - MANOMETER


- Unit of measure selection


Press the **MANOMETER** pane to access the page where you can select the suction pressure unit of measure.

Scale

Allows changing the unit of measure adopted for the suction pressure. To change the unit of measure, press on the arrows until reaching the desired unit of measure. When setting is complete, confirm with **OK**. The main menu will automatically appear.

The operating page will show the pressure with the unit of measure selected.

Note

- *If the unit of measure is changed, the value shown on the scale on the main page will also change.*
- *The values on the various scales are:*
 - e) mmHg: 20 to 200 with increments of 5 units
 - f) cmH₂O: 25 to 265 with increments of 5 units
 - g) mBar: 25 to 265 with increments of 5 units
 - h) Pa x 100 (Pascal x100): 25 to 265 with increments of 5 units

3 – TIMER


- Timer activation and setting

Press the **TIMER** pane to access the timer activation and setting page.

Activating the timer and setting a reminder

To activate this function, press on the switch-shaped button until it turns green and the word **YES** appears. When the function is active, you can change the time by pressing on the relative arrows.

When setting is complete, confirm with **OK**. The main menu will automatically appear. Press **ESC** to return to the operating page where the clock symbol is displayed to show that the timer is active.

Warning

- *The timer has a reminder function only. When the set time has run out, an acoustic signal will be emitted and a message appears on the display to reset it by pressing OK.*
- *When the timer runs out, the functions set on the device remain active (e.g. suction pressure). The device DOES NOT stop performing the functions set.*
- *The countdown of the set time starts from when the data is confirmed.*
- *The timer automatically resets when the device is switched off.*
- *The maximum settable value is 8 hours and the minimum value 10 minutes with increments of 10 minutes.*

4 – ALARMS


- Automatic alarm reset
- Alarm sound volume

Press the **ALARMS** pane to access the page that includes the alarm adjustment functions.

Automatic alarm reset

To deactivate this function, press on the switch-shaped button until it turns red and the word **NO** appears. When setting is complete, confirm with **OK**. The main menu will automatically appear. Press **ESC** to return to the operating page.


- This function enables/disables automatic resetting of the alarms.
- Setting **NO**, all the alarms shown on the main page will be reset only when pressing **OK**.
- Setting **YES**, the alarms will automatically be reset when the fault has been cleared.

Volume


Note

This function allows adjusting the acoustic alarm volume to 10 different increasing levels: 1 = minimum / 10 = maximum The preset volume level is 5.

To activate this function, press on the respective arrows to decrease or increase the alarm volume. When setting is complete, confirm with **OK**. The main menu will automatically appear. Press **ESC** to return to the operating page.

5 – HISTORY


- Events summary
- Events details
- Download events

The **HISTORY** menu contains the data logger function, i.e. you can record the set parameters and the operating events of the device, which can be:

- viewed on the display;
- downloaded to the USB pen drive provided in CSV format and viewed on a PC with Microsoft Excel.


Press the **HISTORY** pane to access the page that includes the events summary, details and download functions.


Note

The parameters recorded in the data logger are deleted if a new therapy is started when switching on the device.

If you unintentionally access the NEW THERAPY menu when switching on the device, you can:

- return to the initial menu by pressing on the **NO** button and then select **CONTINUE THERAPY** by pressing on the button, or

- download the data to a USB pen drive before starting a new therapy (see the paragraph **HISTORY**)


Caution

- *The memory capacity of the data logger allows recording at least 10,000 events; when the memory is full, the oldest data will be overwritten.*
- *To delete the chronology recorded in the data logger, switch off the device with the buttons unlocked, switch it back on and start a new therapy.*

EVENT SUMMARY

Provides the current therapy times displaying 5 records divided into 2 pages:

Therapy start	Provides the date and time the device was switched on and NEW THERAPY mode was selected on the start page.
Total therapy time	Provides the total time of the current therapy from the start of NEW THERAPY .
ON time	Provides the total time of the current therapy in which the device was on without alarms.
OFF time	Provides the total time of the current therapy in which the device was off and/or in alarm (i.e. the therapy interruption time).
ALARM time	Provides the total time of the current therapy in which the device was in alarm (this time is part of the OFF Time).


After viewing the desired pages, press **ESC** to return to the page that includes the events summary, events details and download events functions.

EVENT DETAILS




Provides in detail with date and time all the events and parameters set from the start of the current therapy, including all the alarms that have appeared and been cleared. The display is divided into several pages for easier understanding.

After viewing the desired pages, press **ESC** to return to the page that includes the events summary, events details and download events functions.

Note


In order to correctly read the alarms recorded by the data logger, refer to the table below.

<i>Alarm description</i>	<i>Image on the display</i>	<i>Description in the data logger</i>
INSERT SD CARD! CALL SERVICE		Insert SD card alarm
PUMP FAILURE! CALL SERVICE		Pump malfunction alarm
FULL CANISTER ALARM!		Canister alarm
EXHAUSTED BATTERY! CONNECT THE BATTERY CHARGER		Exhausted battery
PRESSURE ALARM! CHECK CIRCUIT		Pressure alarm
CIRCUIT OCCLUSION ALARM! CHECK PATENCY		Tube alarm
LOW BATTERY! CONNECT THE BATTERY CHARGER		Low battery

EVENT DOWNLOAD

Provides the possibility of downloading the data stored/recoded by the device to a USB pen drive (type: "Eurosets" in blue) provided.


A file in CSV format will be created on the pen drive, which can be viewed and changed using Microsoft Excel. This file contains all the device data (reference, serial number and registration) as well as the data recorded in **EVENT SUMMARY** and **EVENT DETAILS**.

Insert the USB pen drive provided with the device.


Then follow the on-screen instructions:

Remove the pen drive and insert it in a PC that has Microsoft Excel installed.

The pen drive will have a file in CSV format.

The filename is the serial number of the device to which the data was downloaded.

When you open the file, you will find all the data set and recorded by the device divided by type:

The device data are shown at the top: REF.... / SN.... / REG....

The therapy times, i.e. the times viewable in the **EVENT SUMMARY** on the device are shown in the centre.

All the events recorded and viewable in the **EVENTS DETAILS** are shown at the bottom.

	Appunti	Carattere
1	REF	EU3902
2	SN	29008
3	M	1
4		
5	START THE	12 Dec 14 16:41
6	TOTAL TIME	2 dd 23 hh 5 mm
7	TIME ON	0 dd 0 hh 17 mm
8	TIME OFF	2 dd 22 hh 47 mm
9	TIME IN ALARMS	0 hh 0 mm
10		
11	DATE	TIME EVENT
12	12 Dec 14	16:43 Unit OFF
13	12 Dec 14	16:44 Unit ON
14	12 Dec 14	16:44 Unit OFF
15	12 Dec 14	16:44 Unit ON
16	12 Dec 14	16:44 Unit OFF
17	12 Dec 14	16:44 Unit ON
18	12 Dec 14	16:44 Unit OFF
19	12 Dec 14	16:44 Unit ON
20	12 Dec 14	16:44 Unit OFF
21	12 Dec 14	16:44 Unit ON
22	12 Dec 14	16:45 Container alarm
23	12 Dec 14	16:45 Unit OFF
24	12 Dec 14	17:28 Unit ON
		17:29 Unit OFF


ITEM IN EXCEL
THERAPY START
TOTAL THERAPY TIME
ON TIME
OFF TIME
ALARM TIME
DATE
TIME
EVENT
Jan
Feb
Mar
Apr
May
Jun
Jul
Aug
Sep
Oct
Nov
Dec
dd
hh
min

Note


In order to correctly read the alarms recorded by the data logger, refer to the table in the paragraph EVENTS DETAILS.

To delete the chronology recorded in the data logger, switch off the device with the buttons unlocked, switch it back on and start a new therapy.

6 – KEY LOCK

- Touchscreen button locking

This function allows locking the buttons that appear on the screen to prevent voluntarily or involuntarily changing parameters or settings.

Caution


- When switching the device off and on again with the buttons locked, the start page does NOT appear, but it comes on in CONTINUE THERAPY mode, hence with the same settings/parameters as before switching off.
- In addition, in this mode, the data recorded by the data logger are NOT cancelled but recording is continued.

KEY LOCK

When this function is active, all the display buttons are locked, except the **MENU** button on the operating page, the **KEY LOCK** pane and the physical on/off button of the device.

To activate button locking, press the **KEY LOCK** pane to access the relative page where you can enter the following four-digit code using the arrows: 2015. Confirm by pressing **OK**.

Button locking is indicated by a padlock that appears on the operating page.


The buttons remain locked even if the device is switched off and on again.


Caution

- When switching the device off and on again with the buttons locked, the start page does NOT appear, but it comes on in CONTINUE THERAPY mode, hence with the same settings/parameters as before switching off.
- In addition, in this mode, the data recorded by the data logger are NOT cancelled but recording is continued.

To unlock the buttons, repeat the locking sequence (numerical code 2015):


7 - SETTINGS

- Device settings menu

Press the **SETTINGS** pane #7 to access the device settings menu.

7.1	LANGUAGES: device language setting page
7.2	DATE/TIME: device date and time setting page
7.3	DISPLAY: LCD touchscreen display setting page
7.4	ALARM CHECK: functional test of the visual and acoustic alarm system
7.5	RESET: initial setting reset page
7.6	MAINTENANCE: maintenance menu accessible only by the device manufacturer
7.7	UPDATE: device software update page

7.1 - LANGUAGES

- Language selection:
ITALIANO
ENGLISH
DEUTSCH
FRANÇAIS
ESPAÑOL
TÜRKÇE

Pressing on the **LANGUAGES** pane you access a page where you can change the device communication language. To select the language, press on the flag until you reach the desired one, then press **OK** to confirm. The **SETTINGS** menu will automatically appear.

7.2 - DATE / TIME


- Current date/time change

Pressing on the **DATE/TIME** pane you access a page where you can change the device date and time. Use the arrows to change the date and time, then press **OK** to confirm. The **SETTINGS** menu will automatically appear.

Caution

The device date/time appears in the HISTORY, therefore the date and time should be correct for purposes of the procedure.

7.3 - DISPLAY

- Keyboard beep
- Screensaver

Pressing on the **DISPLAY** pane you access a page where you can change the display and touchscreen settings.


Keyboard beep

This function enables/disables sounding of a beep each time a button is pressed.


To activate the **Keyboard beep** function, press on the switch-shaped button until it turns green and the word **YES** appears.

When setting is complete, confirm with **OK**. The **SETTINGS** menu will automatically appear.

Screensaver

This function allows automatically putting the screen to sleep after a predefined time in order to save battery life. The function is active by default with a pre-set switch-off time; you can set a different **Delay** time for shutdown. To disable the **Screensaver** function, press on the switch-shaped button until it turns red and the word **NO** appears.

To reactivate the function, repeat the process in the reverse order.

Then press the arrows to set the desired **Delay** time after which the screen will go to sleep.

When setting is complete, confirm with **OK**. The **SETTINGS** menu will automatically appear.

Note

When the screensaver is active, only the device screen stays off, while the green LED located above the on/off button stays on to indicate that the therapy is continuing normally.

To wake up the screen, press anywhere on it.

7.4 - ALARM CHECK


- Functional test of the visual and acoustic alarm system

This function allows testing the acoustic and visual alarms.

This test must be carried out by healthcare workers the first time the device is switched on and before transferring the device to another patient.

To carry out the **ALARM CHECK**, operate as described in the paragraph "Checking the alarms".

7.5 - RESET

- Initial setting reset

Pressing on the **RESET** pane you access a page where you can reset the initial settings (factory settings) of the device.

Preset or automatically resettable parameters each time the **RESET** function is activated:


- *Timed intermittent operation: not active*
- *Default pressure: 50 mmHg*
- *Suction pressure unit of measure: mmHg*
- *Timer (reminder): not active*
- *Alarm reset: automatic when the fault has been cleared*
- *Button beep: not active*
- *Language: English*
- *Alarm volume: level 5*
- *Screen Saver : active*

Note

Resetting the default parameters does not delete the HISTORY of the data logger!

The parameters recorded in the data logger are deleted only when a new therapy is started when switching on the device.

To activate the **RESET** function, press on the switch-shaped button until it turns green and the word **YES** appears. When setting is complete, confirm with **OK**. The **SETTINGS** menu will automatically appear.


7.6 - MAINTENANCE

THIS FUNCTION IS NOT ACCESSIBLE TO THE OPERATOR BUT ONLY TO THE DEVICE MANUFACTURER.

If you access the **MAINTENANCE** menu by mistake, when the password entry page appears, press **ESC** to return to the **SETTINGS** menu.


7.7 - UPDATE




- Software update



If a new software update becomes available, this function allows updating the device software via USB. Pressing on the **UPDATE** pane you access a page where you can update the software. Insert the USB pen drive and confirm by pressing **OK**. When the update is complete, press **ESC** to return to the **SETTINGS** menu.



Alarms



The WaterLily Vacuum Generator TS device is equipped with the following visual and acoustic alarms:




Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
PUMP FAILURE! CALL SERVICE activated if a functional pump fault is detected	 The alarm page disappears when the problem is solved or after pressing OK, the operating page appears cyclically every 15 sec showing the alarm silenced symbol and the alarm present symbol blinking in red  and for 5 seconds  If the alarm condition has not been cleared 2 minutes after pressing OK, the alarm page reappears permanently 	On (device operational)	Blinking red LED on	On Silenceable for 2 minutes by pressing OK	<p>Switch off the device and request technical service by the manufacturer or an authorised agent.</p> <p> WARNING The PUMP FAILURE! CALL SERVICE alarm stops the device, with consequent blocking of the aspiration.</p> <p> WARNING The absence of vacuum in the dressing might determine build-up of exudates triggering infectious and inflammatory processes. If the absence of vacuum extends for several days, necrotic eschar may form. If suction is interrupted for more than 2 hours, remove the dressing after irrigating it with saline solution.</p>


Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
INSERT SD CARD! CALL SERVICE activated if switching on the device, the absence of the SD memory card is detected.	 <p>The alarm page disappears when the problem is solved or after pressing OK, the operating page appears cyclically every 15 sec showing the alarm silenced symbol and the alarm present symbol blinking in red</p>  <p>and for 5 seconds</p>  <p>If the alarm condition has not been cleared 2 minutes after pressing OK, the alarm page reappears permanently</p> 	On (device operational)	Blinking red LED on	Off	Switch off the device and request technical service by the manufacturer or an authorised agent.

Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
PRESSURE ALARM! CHECK CIRCUIT activated if the set pressure has not been reached within 2 minutes or the pressure remains below the set value for more than 2 minutes or the pressure remains below the set value for three consecutive times in less than 2 minutes	 <p>The alarm page disappears when the problem is solved or after pressing OK, the operating page appears cyclically every 15 sec showing the alarm silenced symbol and the alarm present symbol blinking in yellow</p>  <p>and for 5 seconds</p>  <p>If the alarm condition has not been cleared 2 minutes after pressing OK, the alarm page reappears permanently</p> 	On (device operational)	Blinking Yellow LED on	On Silenceable for 2 minutes by pressing OK	<p>Check that the canister and the vacuum generator are properly connected.</p> <p> WARNING</p> <p>Check that the wound is not haemorrhaging (check for any signs of bleeding in the suction tube and the canister). If so, urgent and appropriate medical action needs to be taken.</p> <p>Check that the wound is not taking in air and/or that the suction tube connections are sealed.</p> <p>The alarm can be silenced for 2 minutes by pressing OK. This alarm is automatically reset only when the problem has been solved (when the set pressure has been reached).</p> <p> WARNING</p> <p>The absence of vacuum in the dressing might determine build-up of exudates triggering infectious and inflammatory processes. If the absence of vacuum extends for several days, necrotic eschar may form. If suction is interrupted for more than 2 hours, remove the dressing after irrigating it with saline solution.</p>
					Check the level of exudates in the collection canister.

Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
POSSIBLE FULL CANISTER ALARM! activated if the fluid collection canister is full	 <p>The alarm page disappears when the problem is solved or after pressing OK, the operating page appears cyclically every 15 sec showing the alarm silenced symbol and the alarm present symbol in red</p>  <p>and for 5 seconds</p>  <p>If the alarm condition has not been cleared 2 minutes after pressing OK, the alarm page reappears permanently</p> 	On (device operational)	Blinking Red LED on	On Silenceable for 2 minutes by pressing OK	<p> WARNING If the canister is full, before replacing it and resuming the therapy, check that the wound is not haemorrhaging (check for any signs of bleeding in the suction tube and the canister). If so, urgent and appropriate medical action needs to be taken.</p> <p>Check that the canister has not tilted/tipped over. If so, move it back into vertical position. If the alarm automatically resets, the canister does not have to be replaced, otherwise, even if the canister is not full, replace it and resume the therapy.</p> <p>The alarm can be silenced for 2 minutes by pressing OK This alarm is automatically reset only when the problem has been solved.</p> <p> WARNING During the alarm silence time, keep the patient and the device monitored to check that the alarm condition has actually been cleared.</p> <p> WARNING The absence of vacuum in the dressing might determine build-up of exudates triggering infectious and inflammatory processes. If the absence of vacuum extends for several days, necrotic eschar may form. If suction is interrupted for more than 2 hours, remove the dressing after irrigating it with saline solution.</p>

Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
<p>CIRCUIT OCCLUSION ALARM! CHECK PATENCY</p> <p>activated if there is an obstruction in the tubes of the suction circuit</p>	 <p>The alarm page disappears when the problem is solved or after pressing OK, a flushing phase starts during which the IN PROCESS FLUSHING... page appears cyclically every 5 sec showing the alarm silenced symbol and the alarm present symbol blinking in yellow</p>  <p>and for 5 seconds</p>  <p>If OK is not pressed within 2 minutes to reset the alarm, the alarm page reappears</p> 	<p>On (device operational)</p>	<p>Blinking Yellow LED on</p>	<p>On Silenceable for 2 minutes by pressing OK</p>	<p>Check that the drain tubes are not kinked and/or squeezed. Remove the cause of kinking/squeezing. Check that any clamps on the drain tubes are open and if not, open them.</p> <p>If none of the above conditions is true and the alarm persists, check that the gauze in contact with the wound has a wrinkled appearance indicating that the wound is being subjected to suction and the tube is hence patent.</p> <p>If the wound is NOT subjected to suction, replace the canister and the drain tubes.</p> <p>The alarm can be silenced for 2 minutes by pressing OK This alarm is automatically reset only when the problem has been solved.</p> <p>WARNING The absence of vacuum in the dressing might determine build-up of exudates triggering infectious and inflammatory processes. If the absence of vacuum extends for several days, necrotic eschars may form. If suction is interrupted for more than 2 hours, remove the dressing after irrigating it with saline solution.</p>

Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
LOW BATTERY! CONNECT THE BATTERY CHARGER activated if the rechargeable internal battery is low	 <p>The alarm page disappears when the power supply/battery charger is connected or after pressing OK, the operating page appears cyclically every 15 sec showing the alarm silenced symbol and the alarm present symbol in yellow</p>  <p>and cyclically every 15 sec.</p>  <p>If the alarm condition has not been cleared 5 minutes after pressing OK, the alarm page reappears permanently</p> 	On (device operational)	Blinking Yellow LED on	On Silenceable for 5 minutes by pressing OK	<p>Power the device directly from the mains using the power supply/battery charger which will at the same time start a new recharging cycle.</p> <p>This alarm can be reset by pressing OK or is automatically reset when the problem has been solved.</p> <p></p> <p>During the alarm silence time, keep the patient and the device monitored to check that the alarm condition has actually been cleared.</p>

Alarm	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/ red LED	Speaker	Action to take
EXHAUSTED BATTERY! CONNECT THE BATTERY CHARGER activated if the rechargeable internal battery is exhausted	 <p>The alarm page disappears when the power supply/battery charger is connected</p>	On (device operational)	Blinking Red LED on	On NOT SILENCEABLE	Power the device directly from the mains using the power supply/battery charger which will at the same time start a new recharging cycle.



Information signals


The WaterLily Vacuum Generator TS device is equipped with the following visual and acoustic information signals:

Information signal	Warnings shown on the WaterLily Vacuum Generator TS display	Green LED	Yellow/red LED	Speaker	Action to take
SEAL TEST activated at the start of therapy until the set suction pressure has been reached		On (device operational)	Off	On a two-tone acoustic signal is sounded when the set suction pressure has been reached.	None.
REMINDER activated when the set reminder time has run out		On (device operational)	Off	On a two-tone acoustic signal is sounded when the time has run out	Reset the reminder signal by pressing OK .

Alarm and information signal priorities in case of overlapping

If more than one alarm occurs at the same time, the device will signal the alarms one at a time with the following priorities:

Type of alarm or signal	INSERT SD CARD! CALL SERVICE	PUMP	POSSIBLE FULL CANISTER ALARM!	BATTERY EXHAUSTED! CONNECT THE BATTERY CHARGER	PRESSURE ALARM! CHECK CIRCUIT
Priority from MAXIMUM to MINIMUM	MAXIMUM PRIORITY! 1	2	3	4	5
Intervention priority from HIGH to LOW	high priority	high priority	high priority	high priority	medium priority
Page	 A.8 INSERT SD CARD! CALL SERVICE OK	 A.7 PUMP FAILURE! CALL SERVICE OK	 A.2 ALARM POSSIBLE FULL CANISTER! OK	 A.0 EXHAUSTED BATTERY! CONNECT THE BATTERY CHARGER OK	 A.1 PRESSURE ALARM! CHECK CIRCUIT OK

Type of alarm or signal	CIRCUIT OCCLUSION ALARM! CHECK PATENCY	LOW BATTERY! CONNECT THE BATTERY CHARGER	SEAL TEST	REMINDER
Priority from MAXIMUM to MINIMUM	6	7	8	MINIMUM PRIORITY! 9
Intervention priority from HIGH to LOW	medium priority	medium priority	information signal	information signal
Page	 A.4 CIRCUIT OCCLUSION ALARM! CHECK PATENCY OK	 A.0 LOW BATTERY! CONNECT THE BATTERY CHARGER OK	 04/08/2015 Tue 10:33:33 80 mmHg LEAK High Low SEAL TEST MENU	 A.3 REMINDER OK

Device classification according to the standards and relative limits of use

Classification according to EN 60601-1 (2006) and EN 60601-1-2 (2007) regarding medical electrical devices and EN ISO 10079-1 (2009) "Medical suction equipment - Part 1: Electrically powered suction equipment - Safety requirements"

- Definition:
 - Mobile/portable device
 - Device for continuous operation

Classification according to EN 60601-1-11 (2010) regarding medical electrical devices used in a domestic environment

- Definition:
 - Transit-operable/body-worn
- Electrical hazards/protection against electric shock
 - Class II medical electrical equipment
 - Medical electrical equipment powered by an internal electric source
- Degree of protection against contacts
 - Device with applied part type BF considered over the entire surface of the device
- Protection against hazardous penetration of solid bodies and water
 - Device with IP22 degree of protection
 - IP40 degree of protection limited solely to the power supply/battery charger provided
- Sterilization methods
 - Device NOT to be sterilized
- Possibility of use in an oxygen-rich environment
 - Device not to be used in an oxygen-rich environment
 - Device not to be used in the presence of flammable substances or agents

Classification according to MDD 93/42/EEC and MDD 2007/47/EC regarding medical devices

- Class IIb device
- The device has passed the tests (where applicable) as required by the following standards:
 - EN 60601-1 (2006): "Medical electrical equipment - Part 1: General requirements for basic safety and essential performance"
 - EN 60601-1-2 (2007): "Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests"
 - EN ISO 10079-1 (2009): "Medical suction equipment - Part 1: Electrically powered suction equipment - Safety requirements"
 - EN 60601-1-8 (2007): "Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems"
 - EN 60601-1-11 (2010): "Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment".

- Limited solely to the power supply/battery charger provided:
 - The power supply/battery charger has passed all the applicable tests according to EN 60601-1 (2006) except those indicated in clause 15.5.3. Referring to clause 15.5.3, the tests passed are those required by CEIEN 60601-1 (1998) indicated in clause 57.9.

Electromagnetic compatibility

<u>ELECTROMAGNETIC EMISSIONS</u>		
WATERLILY VACUUM GENERATOR TS is intended for use in the electromagnetic environment specified below. The customer or the user of WATERLILY VACUUM GENERATOR TS should assure that it is used in such an environment.		
EMISSION TEST	IEC 60601-1-2 TEST LEVEL	ELECTROMAGNETIC ENVIRONMENT GUIDELINES
RF emissions CISPR 11	Group 1	WATERLILY VACUUM GENERATOR TS only uses RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	WATERLILY VACUUM GENERATOR TS is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliant	
Warnings	<p>WATERLILY VACUUM GENERATOR TS is in compliance with IEC EN 60601-1-2 on electromagnetic compatibility for medical electrical equipment. However, as a precaution, the device should not be used in the immediate vicinity of high-power equipment or equipment that, by their nature, emit strong electromagnetic fields. Mobile phones, portable phones or other radio communications equipment used near the device might affect its operation.</p> <p>Should the device necessarily need to be used in the vicinity of other equipment, it should be observed to check that it is operating normally in the configuration in which it is used.</p> <p>Do not use extension cables or adapters on the power cable.</p> <p>The device must be connected to the mains only with the 4m-long power cable provided; the device has been tested with this cable and resulted compliant to the Electromagnetic Compatibility standards. Using cables of different length may cause increased emissions or reduced immunity against radio disturbances.</p>	

ELECTROMAGNETIC IMMUNITY

WATERLILY VACUUM GENERATOR TS is intended for use in the electromagnetic environment specified below. The customer or the user of **WATERLILY VACUUM GENERATOR TS** should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT GUIDELINES
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic pane. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of WATERLILY VACUUM GENERATOR TS requires continued operation during power mains interruptions, it is recommended that WATERLILY VACUUM GENERATOR TS be powered from an uninterruptible power supply (UPS) or buffered by emergency systems (batteries, generators with combustion engines).
Power frequency magnetic field	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.
Note	U_T is the AC mains voltage prior to application of the test level.		

ELECTROMAGNETIC IMMUNITY

WATERLILY VACUUM GENERATOR TS is intended for use in the electromagnetic environment specified below. The customer or the user of **WATERLILY VACUUM GENERATOR TS** should assure that it is used in such an environment.

<p><i>Conducted RF</i> IEC 61000-4-6</p>	<p><i>3V eff</i> <i>150 KHz to 80 MHz</i></p>	<p><i>3V eff</i></p>	<p>Portable and mobile RF communications equipment (such as domestic wireless devices, mobile phones, cordless phones and their base radio stations, walkie-talkies) should be used no closer to any part of WATERLILY VACUUM GENERATOR TS, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1.2\sqrt{P}$ 150 KHz to 80 MHz</p>
<p><i>Radiated RF</i> IEC 61000-4-3</p>	<p><i>3 V/m</i> <i>80 MHz to 2.5 GHz</i></p>	<p><i>3 V/m</i></p>	<p>$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2.4\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p><i>Interference may occur in the vicinity of equipment marked with the following symbol:</i></p> 
<p>Note</p>	<p>These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which WATERLILY VACUUM GENERATOR TS is used exceeds the applicable RF compliance level above, WATERLILY VACUUM GENERATOR TS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating WATERLILY VACUUM GENERATOR TS.</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3 V/m in the frequency range between 150 kHz and 80 MHz.</p>		

**RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF
COMMUNICATIONS EQUIPMENT AND WATERLILY VACUUM GENERATOR TS**

WATERLILY VACUUM GENERATOR TS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the WATERLILY VACUUM GENERATOR TS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and WATERLILY VACUUM GENERATOR TS as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT OF TRANSMITTER (W)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (m)		
	150 KHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

Note	<p>For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be determinated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>
------	--

Warranty conditions

WATERLILY VACUUM GENERATOR TS is guaranteed for **12 months** from the date of purchase on the condition that the device has not been modified and/or tampered with by unauthorised persons.

THE WARRANTY DOES NOT COVER:

- e) Failures and faults caused by improper use of the device, tampering or negligence.
- f) Operations performed for presumed defects or modifications made for convenience.
- g) Parts subject to normal wear.
- h) Transport, delivery and collection.

In accordance with the above mentioned conditions, in the event of malfunctioning, send the device together with the warranty certificate and a copy of the delivery note or the purchase invoice to:

EUROSETS S.r.l.

Strada Statale 12, 143 - 41036 Medolla (MO) ITALY

The device must be shipped carriage free. EUROSETS S.r.l. does not assume any responsibility for damage to the device or its loss during transport.

Explanation of the symbols used on the device labels

	ON/OFF button
	Acoustic alarm silenced
	Warning!
	Medium priority alarm!! Action required within 10 minutes
	High priority alarm!!! Action required within a few minutes
	Obligation!
	USB port for connection of the pen drive provided
	Use a power supply of the type described <i>Type: 7362M/18</i>
	CLASS II device
	APPLIED PART TYPE BF
	Manufacturer
	Year of manufacture
	Code number
	Serial number
	Collect separately as electrical equipment
	Do not use if package is damaged
	Read the operating instructions
	Maximum storage or operating temperature
	Maximum storage humidity
	Prohibition! Do not dispose of the device with the household waste, but return it to the distributor/manufacturer
	Do not wet Do not expose to the elements
	Protected against: - penetration of solid bodies > 12 mm. - dripping water with maximum 15° inclination
	Protected against penetration of solid bodies > 1 mm <u>Not protected</u> against dripping water.

Medolla (MO) - Italy www.eurosets.it

Medium Vacuum
Low Flow
max -200 mmHg

Vacuum Generator TS

Manufacturer:

EUROSETS S.r.l. 41036 MEDOLLA (MO) Italy
Tel.: +39- 0535-660311 - Fax: +39- 0535-51248
e-mail: info@eurosets.com - homepage: www.eurosets.com

WaterLily™ Wound Therapy
is a trademark of Eurosets S.r.l.

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ist eine Handelsmarke von Eurosets S.r.l.